

AED Trainer | XFT-120C |

Introduction

The XFT-120C AED Trainer with Replaceable Language/Scenario Card is designed to effectively train emergency responders in the use of an Automated External Defibrillator (AED). This training unit provides a variety of pre-programmed scenarios to help responders become familiar with an AED and allows them to demonstrate the basic skills necessary to use an AED in an emergency.

XFT-120C AED Trainer provides simulated shock delivery. It has no high-voltage capabilities, and cannot deliver a shock, thereby ensuring safety during training



Functions and Features

- It is completely compliant with AHA new guidelines, and preconfigured with 10 training scenarios that simulate realistic sudden cardiac arrest episodes.
- It includes the DC and AUDIO interfaces, the former can connect the adapter (that's to say, XFT® AED Trainer can be operated by not only batteries but also the adapter), while the later AUDIO OUT can connect the speaker, it's very helpful for instructors during teaching.
- With 2 different connectors (Yellow & Pink), it can be used with manikins or used separately, both of these two connectors can give response to check if the pads connection is firm or poor.
- It can be set bilingual and any language is available, also has Remote Control function.



Currently there are 16 languages, if another languages are needed, that will be available.

Voice language

English (British English & American English)

French (French accent & Canadian accent)

Spanish (Spanish accent & American accent)

German

Finnish

Polish

Greek

Portuguese

Italian

Croatian

Dutch

Hungarian

Hebrew

Cantonese

Arabic

Japanese